


# **The Sequence Analysis Toolkit for the advanced analysis of EthoVision XT data**



Patrick Zimmerman, Wil van Dommelen  
Noldus Information Technology  
Tuesday August 28, Measuring Behavior 2012

## Patrick Zimmerman

- Behavioral Research Consultant at Noldus IT
  - Consultant for behavioral research
  - Documentation specialist
  - Trainer
  
- Ask questions during Tutorial / contact me at [patrick@noldus.nl](mailto:patrick@noldus.nl)


## **Wil van Dommelen**

- Wrote the Sequence Analysis Toolkit
- EthoVision project leader at Noldus IT

### **Sequence Analysis Toolkit is a useful tool to**

- Analyze the sequence of zone visits / behavioral states in EthoVision XT export data files
- Find specific sequences of zone visits / behavioral states
- Get additional information about EthoVision's behavioral states
- Allows export of sequences of zones / behavioral states to The Observer XT

## Applications of the Sequence Analysis Toolkit



# When do I use the Sequence Analysis Toolkit?

What's the  
answer to Life,  
the Universe  
and Everything?


- When you have questions that EthoVision XT cannot answer; Questions that involve EthoVision XT *state variables* (for example, in Zone, Moving/Not moving) that can translate into *sequences*

Radial maze – How many times were all 8 arms consecutively visited once?

Y-maze – How many 3-arm sequences (1-2-3, 2-3-1, 3-1-2, etc.) were made?


Behavioral state – How many times was the animal not moving in a specific zone for at least 30 seconds (“which means it must have been eating”)?

42!



# When do I use the Sequence Analysis Toolkit?

- When you want to export behavioral states from EthoVision XT to The Observer XT



## EthoVision XT

- Per zone, behavioral state –  
Frequency, Total / Mean duration, Latency to first, Latency to last
- Transition – From Zone A to Zone B

## Sequence Analysis Toolkit


- Per zone, behavioral state –  
Frequency, Minimum / Maximum / Average Latency,  
Minimum / Maximum / Average Duration
- Sequence of zones (e.g., arms, quadrants),  
behavioral states (moving/not moving)



## Sequence Analysis Toolkit

- Is a macro in Excel (version 2007 and later)
- Which reads and processes EthoVision XT export files
- To analyze sequences of Zones and Behavioral states
  
- It is not an 'official' Noldus product
- It is supported by Noldus' behavioral research consultants


## 1. Acquire trials



## 1. Acquire trials

				User-defined	User-defined	User-defined	System	User-defined
Label				Subject number	Exposure	Treatment	Acquisition status	Day
Description							The current status of acquisition per arena	
Type				Text	Text	Text		Text
Format								
Predefined Values					Acquisition; +	None; Vehicle;	Unknown; Postpo	Day 1; Day 2; D
Scope				Subject	Subject	Subject	Arena	Subject
Trial	Arena	Subject	No.					
Trial 1	Arena 1	Subject 1	1	333 1	+ 48 hours	None	Acquired	Day 6
Trial 2	Arena 1	Subject 1	2	333 2	+ 48 hours	Vehicle	Acquired	Day 6
Trial 3	Arena 1	Subject 1	3	333 3	+ 48 hours	Scopolamine	Acquired	Day 6
Trial 4	Arena 1	Subject 1	4	333 4	+ 48 hours	None	Acquired	Day 6
Trial 5	Arena 1	Subject 1	5	333 5	+ 48 hours	Vehicle	Acquired	Day 6
Trial 6	Arena 1	Subject 1	6	333 6	+ 48 hours	Scopolamine	Acquired	Day 6
Trial 7	Arena 1	Subject 1	7	333 7	+ 48 hours	None	Acquired	Day 6
Trial 8	Arena 1	Subject 1	8	333 8	+ 48 hours	Vehicle	Acquired	Day 6
Trial 9	Arena 1	Subject 1	9	333 9	+ 48 hours	Scopolamine	Acquired	Day 6
Trial 10	Arena 1	Subject 1	10	333 10	+ 48 hours	None	Acquired	Day 6
Trial 11	Arena 1	Subject 1	11	333 11	+ 48 hours	Vehicle	Acquired	Day 6
Trial 12	Arena 1	Subject 1	12	333 12	+ 48 hours	Scopolamine	Acquired	Day 6
Trial 13	Arena 1	Subject 1	13	333 13	+ 48 hours	None	Acquired	Day 6
Trial 14	Arena 1	Subject 1	14	333 14	+ 48 hours	Vehicle	Acquired	Day 6
Trial 15	Arena 1	Subject 1	15	333 15	+ 48 hours	Scopolamine	Acquired	Day 6

1. Acquire trials
2. Select Zones in an Analysis Profile



1. Acquire trials
2. Select Zones in an Analysis Profile
3. Export the raw data to Excel or Plain text

Raw Data Export

*Data to Export*

☒ Track & dependent variables

☐ Manual scoring log

☐ Trial Control log

☐ Hardware log

*Export Settings*

Destination folder: D:\Temp\SAT Input Browse...

File type: Plain text Delimiter: ;

Start export

## Basic procedure

### Sequence Analysis Toolkit


4. Prepare files to be processed with the Sequence Analysis Toolkit
5. Select files and parameters in the Analysis Profile and set criteria to find specific sequences
6. Carry out the sequence analysis
7. Look at the analysis results

## Sub-sequence *Producer* → *Examiner*

### *Producer*

1. Define a *FromZone* and a *ToZone*
  - Track boundary (first / last zone)
  - A user-specified zone
  - Each zone that occurs in the track
  - A specific sequence or a list of specific zones or sequences
  - FromZone = ToZone
2. Define how to treat multiple encounters of FromZone and ToZone

## Example – Y-maze





### Complete sequence

[**3-3-3-2-1-3-1-2-2-3-1-2-3-1-2-1-3-**]


Let FromZone = **1** and ToZone = **2**

[**3-3-3-2-1-3-1-2-2-3-1-2-3-1-2-1-3-**] ?

[**3-3-3-2-1-3-1-2-2-3-1-2-3-1-2-1-3-**] ?

[**3-3-3-2-1-3-1-2-2-3-1-2-3-1-2-1-3-**] ?

How to treat multiple encounters of FromZone and ToZone?




## Multiple encounters of FromZone and ToZone

[3-3-3-2-1-3-1-2-2-3-1-2-3-1-2-1-3-]

Let FromZone = 1 and ToZone = 2

- **Longest** – First encountered FromZone to last encountered ToZone  
[3-3-3-2-1-3-1-2-2-3-1-2-3-1-2-1-3-]
- **Shortest** – Last encountered FromZone to first encountered ToZone  
[3-3-3-2-1-3-1-2-2-3-1-2-3-1-2-1-3-]  
[3-3-3-2-1-3-1-2-2-3-1-2-3-1-2-1-3-] ... etc.
- **All Sequences** – Each FromZone to each ToZone  
[3-3-3-2-1-3-1-2-2-3-1-2-3-1-2-1-3-]  
[3-3-3-2-1-3-1-2-2-3-1-2-3-1-2-1-3-]  
[3-3-3-2-1-3-1-2-2-3-1-2-3-1-2-1-3-]  
[3-3-3-2-1-3-1-2-2-3-1-2-3-1-2-1-3-]  
[3-3-3-2-1-3-1-2-2-3-1-2-3-1-2-1-3-]  
[3-3-3-2-1-3-1-2-2-3-1-2-3-1-2-1-3-] ... etc.




## Sub-sequence *Producer* → *Examiner*

### *Examiner*


**inspects each sub-sequence given by the *Producer***

- Sub-sequence must exactly match with a specific zone, sequence or list of zones/sequences
- Zones may/may not be revisited in the sub-sequence
- Sub-sequence must **include** a specific zone, sequence or list of zones/sequences
- Sub-sequence must **exclude** a specific zone, sequence or list of zones/sequences
- Subsequence must contain more than, less than or exactly a specific number of zones

**and/or**




## DEMO – Y-maze



**This finds the complete sequence of arm-visits**


Range Criteria:			Define your named lists here:	
From zone:	Track boundary		List 1	List 2
To zone:	Track boundary		any one of	any one of
Focus on:	longest only			
<b>Additional Criteria:</b>				
Exact sequence:				
Pass through ?	visit each zone any nr of times			
Including:				
Excluding:				
Sequence must be:		zones		



[AnalysisProfile.xlsx](#)  
[Results.xlsx](#)

## This finds the alternations in the Y-maze


Range Criteria:			Define your named lists here:	
From zone:	Each zone		List 1	List 2
To zone:	Each zone		any one of	any one of
Focus on:	shortest only		1-2-3	
			1-3-2	
			2-1-3	
			2-3-1	
			3-2-1	
			3-1-2	
<b>Additional Criteria:</b>				
Exact sequence:	List 1			
Pass through ?	visit each zone any nr of times			
Including:				
Excluding:				
Sequence must be:		zones		



[AnalysisProfile.xlsx](#)  
[Results.xlsx](#)


## This finds all the right turns in the Y-maze

Range Criteria:			Define your named lists here:		
From zone:	Each zone			List 1	List 2
To zone:	Each zone			any one of	any one of
Focus on:	shortest only			2-1	
				1-3	
				3-2	
<b>Additional Criteria:</b>					
Exact sequence:	List 1				
Pass through ?	visit each zone any nr of times				
Including:					
Excluding:					
Sequence must be:			zones		



[AnalysisProfile.xlsx](#)  
[Results.xlsx](#)

## DEMO – Radial maze



**Excited,  
are you?!**





## Radial maze


- Short-term memory – How do animals search the for baited arms?  
Do they revisit arms?
- Treatments – How do different treatments (drugs, housing, stress) affect performance in the maze and therefore memory?

**This finds the sequence in which all 8 arms have been exactly once without re-visits**

Range Criteria:				Define your named lists here:	
From zone:	Each zone			List 1	List 2
To zone:	Each zone			any one of	any one of
Focus on:	longest only				
<b>Additional Criteria:</b>					
Exact sequence:					
Pass through ?	visit each zone only once				
Including:					
Excluding:					
Sequence must be:	exactly	8	zones		

[AnalysisProfile.xlsx](#)


[Results.xlsx](#)



**How many times did I visit each arm exactly once in a specific order?**

Range Criteria:		Define your named lists here:	
From zone:	Track boundary	List 1	List 2
To zone:	Track boundary	any one of	any one of
Focus on:	longest only	1-2-3-4	
<b>Additional Criteria:</b>			
Exact sequence:	List 1		
Pass through ?	visit each zone any nr of times		
Including:			
Excluding:			
Sequence must be:			zones


[AnalysisProfile.xlsx](#)  
[Results.xlsx](#)



## How many times did I re-visit an arm?


Range Criteria:			Define your named lists here:	
From zone:	From zone equals To zone		List 1	List 2
To zone:	From zone equals To zone		any one of	any one of
Focus on:	shortest only			
<b>Additional Criteria:</b>				
Exact sequence:				
Pass through ?	visit each zone any nr of times			
Including:				
Excluding:				
Sequence must be:	exactly	2	zones	

[AnalysisProfile.xlsx](#)  
[Results.xlsx](#)



## How many times did I re-visit a specific arm?

Range Criteria:			Define your named lists here:	
From zone:	From zone equals To zone		List 1	List 2
To zone:	From zone equals To zone		any one of	any one of
Focus on:	shortest only		1-1	
			3-3	
			5-5	
			7-7	
<b>Additional Criteria:</b>				
Exact sequence:	List 1			
Pass through ?	visit each zone any nr of times			
Including:				
Excluding:				
Sequence must be:		zones		



## Dreamy 'state'

- 'Dreamy' defined as the animal not moving for at least 10 seconds as defined in EthoVision XT



## DEMO – Behavioral state 'Dreamy'


**Surely, you  
must be fast  
asleep by now?**



# Example – Exporting Zones / Behavioral states to The Observer XT

## EthoVision XT

1. Acquire trials
2. Select Movement states in an Analysis Profile
  - Set Threshold
  - Select all states
3. Export the raw data to Excel or Plain text



**Movement**

Movement Body Points Statistics

Calculates the duration for which the selected body point was changing location.

Filter

Averaging interval: 1 samples

Threshold

Start velocity: 10.00 cm/s

Stop velocity: 5.00 cm/s

Calculate statistics for

☒ Moving

☒ Not moving

OK Cancel





## How many times was the animal 'dreamy' for at least 10 seconds?

Range Criteria:			Define your named lists here:	
From zone:	D		List 1	List 2
To zone:	A		any one of	any one of
Focus on:	shortest only			
Additional Criteria:				
Exact sequence:				
Pass through ?	visit each zone any nr of times			
Including:				
Excluding:				
Sequence must be:		zones		




**AnalysisProfile.xlsx**


## Results.xlsx

## Sequence Analysis Toolkit

1. Click the *Export Event Logs* button
2. Import the export file (\*.txt) into The Observer XT



## Y-maze - Alternations



## **Sequence Analysis Toolkit is a useful tool to**

- Analyze the sequence of zone visits / behavioral states
- Find specific sequences of zone visits / behavioral states
- Get additional information about EthoVision's behavioral states
- Allows export of sequences of zones / behavioral states to The Observer XT



**SO LONG  
&  
THANKS FOR ALL THE FISH**



**Glad this is  
over...**

**Questions?**

